

DERMA

THE RIGHT TOOL FOR MICRODERMABRASION

A blackhead remover that unclog pores



- Charging Voltage > 5V/1A
- Rated Output Power > 0.75W
- Battery Voltage > 3.7V
- Rated Current > 200mA

● DERMA Glacier White

EXFOLIATING AND BLACKHEAD REMOVING WITH PROFESSIONALISM

Proceed microdermabrasion with the DERMA reduces or eliminates the wrinkles and fine lines, "crow's feet," age spots, light acne scars, and undesired pigmentation that may appear on the skin over time.

DERMA peels and removes dead skin cells on the top layer and vacuuming/suction to remove the particles along with dirt and blackheads. Skin is left smoother and clearer.



NON-CHEMICAL. NON-INVASIVE

The DERMA microdermabrasion is a non-chemical, non-invasive procedure that use four different probes to remove the outermost layers of dry, dead skin cell and reveal younger, healthier-looking skin.

4 LEVELS OF INTENSITY

Soft, normal, strong, super. Pick an intensity level that suits your skin.

Note: soft level is recommended for first-time users.

4 TYPES OF PROBE. DIFFERENT PURPOSES



Microcrystalline Probe
The microcrystalline probe emits millions of microscopic crystals onto the skin with high absorption to remove dead skin cells. Gentle on skin compared to traditional demabrasion facials.



Large Circular Probe
The large circular probe using strong suction to massage, lift and stimulate microcirculation and blood flow to tighten skin.



Small Circular Probe
The small circular probe can be aimed at blackheads and pimples to help unclog pores and reduce the look of blemishes. Suitable for sensitive skin.



Oval Probe
The oval probe with micro pulling effect that smooths areas with fine lines including the corners of eyes and mouth. It deeply penetrates the skin, stimulating circulation to prevent future wrinkles from forming.

Wirelessly charge the DERMA.



Get a clearer look at your skin with the build-in LED lighting system