

Face and neck feel lifted almost instantly


Stimulate blood flow, skin looks more radiant

More voluminous cheek, and chiseled jawline

Lines and wrinkles become less visible Promote the growth of collagen, skin feel firmer

## A COMBINATION OF MICRO-CURRENT AND HI-FREQUENCY SONIC TECHNOLOGY

The KAXA-LIFT is an ergonomically shaped sonic face massager that helps to refresh your face by targeting the most problematic areas, including your forehead, eye area, brows, liplines, cheeks and chin. All you need is 6 minutes a day to achieve a more radiant, younger looking appearance.





### MICROCURRENT

The MARASIL KAXA-LIFT utilizes microcurrent technology for facial stimulation. The micro-current sends soft, gentle waves through the skin, tissues and down to the facial muscles. It stimulates ATP production which drives the creation of collagen, and give you an improved, lifted appearances.

### **HI-FREQUENCY SONIC**

The Hi-Frequency Sonic pulsation helps to stimulate blood circulation, and boosts the absorption of your favorite daily anti aging skincare.

### ERGONOMICALLY SHAPED TO TARGET THE MOST PROBLEMATIC AREAS



### TWO SPECIALIZED MODES

#### HEATED

- · Facial lifting.
- · Boost the absorption of your daily anti aging skincare.
- · Lines and willides reduc
- · Facial muscle massage.
- · Promote collagen production
- · Stimulate blood circulation.





### NATURAL

· Improve skin tone.
· Reduce Edema and Eye bags.
· Reduce pore size.

Restore skin elasticity and tenderness
fassage acupuncture points to release stress
Promote collagen production
Stimulate blood circulation

### Two specialized mode

Normal and Heated. Normal for tightening and firm skin. Heated for massaging and promote skin absorption

### Microcurrent

Promote the production of collagen

### **Hi-frequency sonic**

Better massaging, better absorption



# **Ergonomically Design From Japan**

Improve your skin at every angle

### **Travel-friendly**

compact and portable, facial lifting whereverand whenever

### Intelligent

Touch button to activate the device; 1 minute timer remind you to target next area; automatically shut off when left on for 2 minutes