

# PLASMA PEN

## 2 IN 1 ANTI-AGING FIBROBLAST



Charging Parameters > DC5V / 500mA  
Product Power > 3W  
Charging Time > 2.5H  
Rated Current > 800mA

# INNOVATIVE TECHNOLOGY FOR HOME USE



Collagen  
Reproduction



Deep  
Absorption



Anti  
Inflammation



Skin  
Resurface



Tightening and  
Face Lifting



Safe and  
Portable

Plasma head

①

30secs per Section


②

10mins everyday

on/off




## PLASMA FIBROBLASTING TECHNOLOGY




## OXYGEN AND NITROGEN OZONE RELEASE

Help Acne Recovery



Ozone sterilization



| Kill Bacteria



| Improve Breakout

