

MICROCURRENT PHOTON THERAPY BEAUTY DEVICE



Rated Power	>	5w
Rated Voltage	>	DC5V
Battery Capacity	>	920mAH
Charging Time	>	2.5hrs
Color	>	White

✉ hello@bolva-co.com

🌐 www.bolva-co.com

EMS MICROCURRENT MODE

Stimulate Subcutaneous Cell Reproduction, Improvement Skin Texture through Effective Lifting and Tightening



Microcurrent Stimulation



ION Penetration and Derivation



Enhance Absorption

LED LIGHT THERAPY



Blue Light Mode

Wavelength: 460-470nm


1. Anti-Inflammation
2. Acne Clearance
3. Facilitate Healing



Yellow Light Mode

Wavelength: 585-595nm

1. Reduce Dark Spots
2. Even Skin Tone
3. Improve Brightness



Red Light Mode

Wavelength: 620-630nm

1. Reproduction
2. Erase Finelines
3. Rejuvenation



Multi-Light Mode

Wavelength: 620-630nm

1. All Skin Types
2. All-in-1 Functions
3. Solves All Problems



TOUCH MASSAGING PROBES



CONTROL PANEL BUTTON



USB CHARGING PLUG